

Automatic classification of intensity-vs-time curves in breast DCE-MRI by K-means clustering and Dynamic Time Warping curve matching.


García-Castro F¹, Alberich-Bayarri A^{1,2}, García-Juan D¹, Martí-Bonmatí L²

¹ QUIBIM S.L., Valencia, Spain

² GIBI 230, La Fe Polytechnics and University Hospital, Valencia, Spain



HOSPITAL LA FE DE VALENCIA
UNIVERSITARI I POLITÈCNIC





- Introduction
- Purpose
- Materials and methods
 - K-means
 - Dynamic Time Warping
- Results
- Conclusion






Univeristy Of Virginia. School of Medicine.



Introduction

Innovación y experiencia al servicio del paciente







- Develop an automatic lesion detection algorithm for breast DCE-MRI.
- Imaging biomarker extraction: volume and eccentricity.
- Automatic generation of a quantitative structured report.

- K-means
- Iterative clustering algorithm.
- Reduce the set of intensity-vs-time curves.


$$\mathcal{L} = \sum_{j=1}^k \sum_{i=1}^n \|x_i^{(j)} - \mu_j\|^2$$





- K-means pre-clustering



K-means
→




- Dynamic Time Warping



$$X = (x_1, \dots, x_N), x_i \in \Phi$$


$$Y = (y_1, \dots, y_M), y_j \in \Phi$$




Optimum alignment path P^*

Minimize: $d: \Phi \times \Phi \rightarrow \mathbb{R} \geq 0$

- DTW voxelwise classification



DTW
→





Report

- Patient data
- 3D Views
- Volume
- Eccentricity
- Average curve




Breast DCE-MR Curve Classifier

Client	Modality	Patient ID	Timepoint	Region	Study date	Study received
	DCE-MRI					



Lesion features



Lesion ID	L1	L2	L3	L4
Volume (cm ³)	2.5	0.25	0.13	0.09
Eccentricity	0.64	0.99	0.56	0.99



©2009-2013 quibim. All rights reserved. This report is the property of quibim. It is not to be distributed, copied, or used in any way without the express written permission of quibim. quibim is a trademark of quibim. All other trademarks are the property of their respective owners.





Case 1





Case 2



- Dynamic Time Warping is a very convenient algorithm for curve matching tasks.
- Automatic detection tools can improve time efficiency in radiology workflows.
- Quantitative structured reporting offers a needed help for radiologists in their diagnostics.

Luis Martí Bonmatí – MD, PhD. GIBI PI and QUIBIM
Founder

Ángel Alberich-Bayarri – PhD. GIBI Director and QUIBIM
CEO



QUIBIM Staff


Fabio García Castro - M.Sc

Rafa Hernández Navarro - B.Sc

David García - M.Sc

Encarna Sánchez - M.Sc

Katherine Wilisch R. - M.Sc





GIBI²³⁰ Staff

Enrique Ruiz Martínez – M.Sc

Amadeo Ten Esteve – M.Sc

Ana Penadés - Adm.



Internship Students

Belén Fos Guarinos

Alfredo Torregrosa Lloret

Carlos Moya Claramunt

Ana Jiménez Pastor

Irene Mayorga Ruiz

